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## POINTS ESSENTIELS

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# vétolabo

PARTENAIRE DE VOS RECHERCHES

## **SM 710 Series Universal Surgical Microscope User Manual V 1.0**

---

1. Precautions .....	2
2. Component Names .....	3
3. Component Installation .....	4
4. Operation .....	6
5. Technical Specifications.....	9
6. Troubleshooting.....	10
7. After-Sales Service.....	11

---

# 1. Precautions

## 1-1 Operation


- (1) Do not expose components directly to sunlight. Place them in a dry, clean environment, avoiding high temperatures and severe vibrations.
- (2) The microscope is a precision instrument. Handle it with care and avoid impacts and collisions during transportation.
- (3) To avoid affecting image clarity, prevent dirt or fingerprints from remaining on the lens surface.
- (4) Do not rotate the left and right zoom handwheels in the opposite direction, as this may cause malfunctions.
- (5) Be careful when installing or removing lenses to prevent them from falling.

## 1-2 Maintenance and care

- (1) Keep all lenses clean. If there is fine dust, use a blower or gently wipe with cotton gauze. If there are oil stains or fingerprints, gently wipe with cotton gauze dipped in a 3:7 mixture of ethanol and ether.
- (2) Do not use organic solvents to clean other surfaces of the microscope, especially plastic surfaces. Use neutral detergents for cleaning.
- (3) Do not disassemble the microscope yourself, as this may affect its performance.
- (4) When not in use, cover the microscope with a dust cover to prevent contact with dust. Store it in a moisture-free place to avoid rust or mold.
- (5) To maintain the microscope's performance, regular inspections are recommended (details can be discussed with your nearest agent).

## 2. Component Names

### 2-1 Component Overview

### 3-3 Installing the CCD Camera and Monitor (Optional)


**Special Note:** This step applies only to trinocular microscopes, which have a CCD Camera Interface on the microscope body, as shown below;

Trinocular microscope

---


Trinocular microscopes with a CCD camera interface can be equipped with a camera and/or monitor, allowing real-time viewing, image capture, and video recording, as shown below.

## 4 Operation

### 4- 1 Using the Lifting and Rotating Functions


Loosen the ④ Locking Handle to allow the microscope stand (body) to move up and down (coarse adjustment) and rotate 360° horizontally, as shown below;


---

Rotate the ③ Focusing Fand-wheel to move the microscope stand (body) up and down (fine adjustment), adjusting the working distance between the objective lens and the observed object for optimal clarity. As shown below;

[Special Note: (1) To adjust the tightness of the focusing mechanism, hold one handwheel and rotate the other. Tightness depends on the rotation direction: clockwise tightens, counterclockwise loosens. (2) Adjusting the tightness of the focusing mechanism prevents the microscope body from sliding down during observation and makes focusing more comfortable.]


Loosen the ⑨ Body Locking Handle to allow the microscope body to rotate 360° , as shown below;

---

## 4- 2 Diopter Adjustment and Focusing

Step 1: Rotate the zoom handwheel to maximum magnification.

Step 2: Rotate the right diopter adjustment ring to the 0 mark.

Step 3: Observe through the right eyepiece. If the image is unclear, rotate the focusing handwheel to clarify the specimen image.

Step 4: Rotate the zoom handwheel to minimum magnification.


Step 5: Observe through the right eyepiece. If the image is unclear, rotate the right diopter adjustment ring ② to clarify the specimen.

Step 6: Rotate the zoom handwheel back to maximum magnification. Observe through the right eyepiece. If the image is unclear, repeat steps 3 to 5 for more precise diopter adjustment.

Step 7: Rotate the zoom handwheel to minimum magnification. Observe through the left eyepiece. If the image is unclear, rotate the left diopter adjustment ring to clarify the image.

## 4-3 Interpupillary Distance Adjustment

Hold the left and right prism housings (as shown below) and rotate them in the direction of the arrows until comfortable for binocular observation.

#### 4- 4 Using the Eyecups

1. For those without glasses, hold the diopter adjustment ring to prevent rotation. Rotate the eyepiece to ensure the eyecup fits well with the observer's eye.
2. For those with glasses, remove the eyecup before observation.

#### 4- 5 Working Area

The area inside the base (below the microscope body) is the working area. You can place the observed object directly here or use a surgical table, as shown below;

# 5 Technical Specifications

Eyepiece	Standard Configuration		Additional Objective Lens					
			0.5X		1.5X		2X	
	Working Distance 100mm		Working Distance 165mm		Working Distance 45mm		Working Distance 30mm	
Magnification	Field of View	Magnification	Field of View	Magnification	Field of View	Magnification	Field of View	
10X/20	7X	28.6	3.5X	57.1	10.5X	19	14X	14.3
	45X	4.4	22.5X	8.9	67.5X	3	90X	2.2
15X/15	10.5X	21.4	5.25X	42.8	15.75X	14.3	21X	10.7
	67.5X	3.3	33.75X	6.7	101.25X	2.2	135X	1.7
20X/10	14X	14.3	7X	28.6	21X	9.5	28X	7.1
	90X	2.2	45X	4.4	135X	1.5	180X	1.1

Notes:

1. The working distance is fixed and does not change with magnification.
2. When using an auxiliary objective lens,  
the total magnification = objective magnification × eyepiece magnification × auxiliary objective magnification.

Eyepiece Field of View

Field of View (mm) =

\_\_\_\_\_

Objective Magnification × Auxiliary Objective Magnification

# 6. Troubleshooting

If the microscope's performance is not fully utilized due to unfamiliarity, the following table provides some solutions.

Issue	Cause	Solution
1. Double images not aligned	Incorrect interpupillary distance adjustment	Adjust the interpupillary distance
	Incorrect diopter adjustment	Re-adjust the diopter
	Different magnifications in left and right eyepieces	Install identical eyepieces
2. Dirt in the field of view	Dirt on the specimen	Clean the specimen
	Dirt on the eyepiece surface	Clean the eyepiece

3. Blurry image	Dirt on the objective lens surface	Clean the objective lens
4. Blurry image during zoom	Incorrect diopter adjustment	Re-adjust the diopter
	Incorrect focusing	Re-focus
5. Focusing handwheel not smooth	Focusing handwheel too tight	Loosen appropriately
6. Microscope body slides down during observation, causing blurry images	Focusing handwheel too loose	Tighten appropriately
7. Cut-off in the field of view during observation or photography	Pull rod not fully extended	Extend the pull rod fully
8. Blurry image on the video monitor during zoom	Incorrect camera focus adjustment	Re-adjust the camera focus using the focus ring on the camera adapter
9. Eye strain	Incorrect diopter adjustment	Adjust the diopter correctly
	Unsuitable lighting brightness	Adjust the light intensity knob