


*Surgical Collection*


# Feeding & Dosing Needles

- Handcrafted, polished stainless steel with smooth olive tip
- Uses range from rectal injections to flushing ear canals
- Universally compatible with industry standard Luer mount
- Autoclavable (wrapped at 121°C for 30 minutes)
- Available in sets of two



Our Feeding and Dosing Needles are intelligently designed for surgical precision.

## Diagram



## Sizes

CURVED	STRAIGHT	LENGTH	GAUGE	TIP SIZE
FNC2211/2	FNS2211/2	1.5" 3.8 cm	22	2.4 mm
FNC2011/2	FNS2011/2	1.5" 3.8 cm	20	2.4 mm
FNC182	FNS182	2" 5 cm	18	2.4 mm
FNC184		4" 10 cm	18	2.4 mm
FNC163	FNS163	3" 7.6 cm	16	3 mm
FNC164	FNS164	4" 10 cm	16	3 mm
FNC143	FNS143	3" 7.6 cm	14	4 mm
FNC144	FNS144	4" 10 cm	14	4 mm
FNC124		4" 10 cm	12	6.4 mm
FNC106	FNS106	6" 15.25 cm	10	6.4 mm
FNC104	FNS104	4" 10 cm	10	6.4 mm
FNCA*	FNSA**			

\* One of FNC182, FNC143, FNC106  
\*\* One of FNS182, FNS143, FNS106